

## HYBRID VASCULAR GRAFTS

Company Name	Product Name	US FDA Indicated Use	Materials Used	Nitinol Reinforced Section		Vascular Graft Section		Introducer Size (F)	Recommended Guidewire (inch)
				Diameter (mm)	Length (cm)	Diameter (mm)	Length (cm)		
Gore & Associates	Gore Hybrid Vascular Graft	Intended for use as vascular prosthesis for replacement or bypass of diseased vessels in patients with occlusive or aneurysmal diseases, in trauma patients requiring vascular replacement, for dialysis access, or for other vascular procedures	Nitinol/ ePTFE	6, 7, 8, 9	5	6	50	N/A for open surgical; 11–14 peelaway for over-the-wire	N/A for open surgical; 0.035 for over-the-wire

ePTFE, expanded polytetrafluoroethylene.