

STANDARD SHEATHS

Company Name	Product Name	Size (F)	Length (cm)	Guidewire Diameter (inch)	Curved or Straight	Radiopaque Tip
Argon Medical Devices, Inc.	Access	4-9	5, 7.5, 11	0.021, 0.032, 0.035, 0.038	Straight	No
	One Piece		11			
	Tuohy-Borst					
	Detachable Valve					
Arrow International, a division of Teleflex	Cath Lab/Interventional Sheaths (adult)	4-9	11	0.035	Straight	Yes
	Cath Lab/Interventional Sheaths (pediatric)	4-6	5, 7.5	0.021		
	Arrow-Flex PSI (With Arrow Raulerson Syringe and Cath-Gard Catheter Contamination Guards; With and Without ArrowGard)	6, 7, 8, 8.5, 9	7, 10	0.025-0.035		
B. Braun Interventional Systems	Intradyn Arterial Introducer Sheath and Elite HV Hemostasis Introducers	3-14	7, 11, 23	0.018, 0.021, 0.035, 0.038	Straight	With or without
Boston Scientific Corporation	Super Sheath Introducers	4, 5, 6, 7, 8, 9	11, 25	With wire: 0.035; without wire: 0.038 compatible	Straight	With or without
	Super Sheath XL Introducers	10, 11, 12, 14				
Cook Medical	Performer: Small Valve	4-8	5, 5.5, 7, 13, 25, 30, 35, 45	0.035, 0.038	Straight	With or without
	Performer: Large Valve	9-18	13, 30, 35, 40, 70, 80	0.038		
	Performer: Radial	4-6	5.5, 7, 13, 23	0.018		No
Cordis Corporation	Brite Tip Interventional Sheath	4-11	5.5, 11, 23, 35, 45, 55, 90	0.035	Straight	Yes
	Avanti Sheath Introducer		5.5, 11, 23			No
Galt Medical Corporation	Micro-Access Elite HV Hemostasis Valve Introducer	4-7	5.5, 12	0.018	Straight	With or without
	ReDial 4 cm	5-7	4	0.018, 0.038		
Gore & Associates	Gore DrySeal Sheath	12, 14, 16, 18, 20, 22, 24, 26	28	0.035	Straight	Yes
Medtronic, Inc.	Input TS	5, 6, 7, 8, 9, 10, 11	7, 11, 23	0.038	Straight	Yes
	Input PS	5, 6, 7, 8, 9	11, 23			
Merit Medical Systems, Inc.	Prelude Sheath Introducers	4, 5, 6, 7, 8	11, 23	0.018, 0.025, 0.035, 0.038	Straight	With or without
	Prelude ACT Sheath Introducers	4.5, 5.5, 6.5, 7.5, 8.5	11	0.035, 0.038		Without

STANDARD SHEATHS (CONTINUED)

Company Name	Product Name	Size (F)	Length (cm)	Guidewire Diameter (inch)	Curved or Straight	Radiopaque Tip		
Oscor Inc.	Adelante (Peel Away)	6, 7, 7.5, 8, 8.5, 9, 9.5, 10, 10.5, 11, 12	16, 23	0.035, 0.038	Straight	No		
	Adelante Sigma (Hemostic Valve)	4, 5, 6, 7, 8, 9, 10, 11	10, 23					
	Adelante Sigma Sol (Hemostic Valve)						Yes	
	Adelante Sigma P (Hemostic Valve)	4, 5, 6, 7, 8	5, 10					
St. Jude Medical, Inc.	Ultimum	4–12, 14, 16, 18, 20, 22	5, 8.5, 12, 23, 25, 30, 35	0.035, 0.038	Straight	Yes, EV line		
	Engage	4, 5, 6, 7, 8, 9	12, 25			No		
Terumo Interventional Systems	Pinnacle Introducer Sheath for Peripheral Applications	4, 5, 6, 7, 8, 9, 10, 11	10, 25 (25 does not come with mini guidewire or in 4 F)	Comes with or without 0.035 guidewire	Straight	No		
	Pinnacle Introducer Sheath for Coronary Applications						Comes with or without 0.038 guidewire; the guidewire for the 4 F is 0.035	
	Pinnacle TIF Tip Introducer Sheath for Peripheral Applications	5, 6, 7, 8	10 with 0.035-inch mini guidewire	Mini guidewire included				
	Pinnacle TIF Tip Introducer Sheath for Coronary Applications		10 with 0.038-inch mini guidewire					
	Pinnacle R/O II Radiopaque Marker Introducer Sheath	4, 5, 6, 7, 8, 9	6, 10, 25	Mini guidewire available with 6 and 10 cm; maximum guidewire is 0.038			Yes	
	Pinnacle R/O II HiFlo	6, 7, 8	4	Comes with or without 0.035 guidewire				
	Glidesheath (indicated for radial access)	4, 5, 6	10, 16, 25	Comes with 0.021 guidewire			Spring/straight	No
	GlideAccess Systems (Micro Accessing Systems)	4, 5	10	Comes with 0.018 guidewire; echogenic or nonechogenic needle			Stainless steel or platinum-tipped nitinol wire	
	Obturator	5, 6, 7, 8, 9	10, 25	Not for use with wire			Fits Pinnacle 5–9-F sheaths; 10- or 25-cm lengths	