

MANIFOLDS

Company Name	Product Name	Number of Ports	On/Off Handles	Port Orientation	Body Style	Pressure	Integrated Transducer	Comments
Argon Medical Systems, Inc.	Manifold	1, 2, 3, 4	On, off	Left, right	Fingergrip, block	500	No	Not provided
Merit Medical Systems, Inc.	Merit Manifolds	2, 3, 4, 5	On, off	Left, right	Block, half body	200, 500	Yes	Not provided
	TRAM (transducer and manifold)	3		Right	Half body, wide			Take pressure readings at heart level
	TRAM-P (transducer and manifold with pressure line)			Take pressure readings at any level during procedure				
Navilyst Medical	NAMIC Compensator Manifold	1, 2, 3, 4	On, off	Left, right	Fingergrip, wide	Medium	Yes	Take pressure readings at any level during procedure
	NAMIC Perceptor Manifold					Medium, high		Take pressure readings at heart level
	NAMIC Manifold	2, 3, 4			Fingergrip, wide, block	No	Not provided	
	NAMIC Traditional/Monitoring Manifold	Total ports: 4, 5, 6, 7; side ports: 2, 3, 4, 5			Off	Block	Medium	All female ports