

STOPCOCKS

Company Name	Product Name	Material	Body Style	Lumen Inner Diameter (inch)	Pressure Ratings (PSI)	Fittings
Argon Medical Devices, Inc.	Stopcocks	Body: clear polycarbonate; handle: acetal	1-, 3-, or 4-way	0.056, 0.093	15 (standard), 1,050 (high)	FLL, rotating collar, MLL, RA, male slip tip
B.Braun Interventional Systems Inc.	High-Pressure Stopcocks	Polycarbonate	1- or 3-way	N/A	1,050	Snap-fit rotating adapter
Merit Medical Systems, Inc.	Merit Marquis Stopcock	Body: clear polycarbonate; handle: acetal	1-, 3-, or 4-way	N/A	50, 200, 500, 1,050	Slip luer, fixed male luer, swivel nut, rotating male adapter
Navilyst Medical	NAMIC Stopcocks	Body: clear polycarbonate; handle: acetal	1-, 3-, or 4-way	0.075	200 (low), 400 (medium), 1,050 (high)	FLL, rotating collar, MLL, RA, male slip tip
	Extension Sets With Stopcocks	Tubing style: pediatric, adult, or IV	Tubing lengths (inch): 4.5, 6, 9, 12, 20, or 33	N/A	N/A	MLL, male slip tip

FLL, female luer lock; MLL, male luer lock.